Composition

Penghawar djambi, eugenol, sodium lauryl sulphate, calcium carbonate, excipients.

Properties

ALVEOXYL is a paste for alveolar surgical dressing after tooth extractions, with Penghawar fibers. The paste promotes haemostasis by compression and protects from superinfection (alveolar osteitis) by its barrier effect.

When introduced into the socket, ALVEOXYL easily adheres to the walls due to its fibrous consistency. With the assistance of the patient's tongue movements, it gradually self eliminates without an intervention by the practitioner as the normal healing of the wound takes place.

ALVEOXYL is very easy to apply and requires no special attention from the practitioner other than observation of the healing process.

Indications

ALVEOXYL is a paste used as dressing in case of dry socket or post-extraction dressing following a difficult or traumatic extraction in patients with history of dry sockets.

Contra-indications

ALVEOXYL is not to be used on patients having a history of allergy to one of the components and on deciduous teeth (children under 12 years of age).

Routine use of the product as a post-extraction dressing is contraindicated.

Instructions for use

Take a small pellet (about 0.20 g) of paste and place it gently into the prepared dental socket.

Do not suture.

Warnings and precautions of use

Patients should be told not to wash their mouth vigorously within the 24 hours following the avulsion.

Do not treat more than 5 teeth in one session.

Xi – irritant

- May cause sensitization by skin contact.
- Irritating to skin
- Avoid contact with skin
- Wear suitable gloves
- This material and its container must be disposed of as hazardous waste.

The cap should always be replaced immediately after use.

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