Cresopate

Paste for dental use

Self-hardening, ready-mixed root canal sealer

Formula

Parachlorophenol	7.36 g
Camphor	11.75 g
Monohydrated zinc sulphate	10 g
Excipients: Vinyl chloride-vinyl acetate copolymer, Zinc oxide, Barium sulphate, Glycerin triacetate q.s.ad.	100 g

Properties

CRESOPATE, root canal filling paste contains two antiseptics :
- Parachlorophenol which combines antiseptic properties and the ability to increase the permeability of the dentinal tubules.
- Zinc sulphate.

It contains no formaldehyde, or derivatives of this product.

The radiopaque excipient is self-hardening. It sets in the root canal, slightly increasing in volume and produces a chalk-like mass, which makes it the ideal choice for sealing points into the canal.

Indications

Root-canal filling after treatment of infected canals.
Antisepsis of root-canals and dentinal tubules.

Method of use

Dental use only.

Dental product for professional use only.

After carefully filling and reaming the canal, introduce the paste either in conjunction with a point, or using a rotary paste filler. Do not extrude the paste beyond the apex.

If root-canal filling is to be completed in the same session, quickly mix one drop of distilled water into the paste before introducing the paste in the canal. This significantly reduces the setting-time.

Important note

CRESOPATE hardens in the presence of moisture, therefore it is important, especially when filling multi-rooted teeth, to **thoroughly dry the pulp chamber** before introducing the paste in the canal.

Use dry instruments for handling the paste.

Close the jar tightly after each use.

Presentation

Jar containing 7.5 g of paste.

Storage

Do not use after the expiry date mentioned on the outer package.

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